

Title:	BCR Height	ID:	
		0012	
Date in:	Response:	Model:	Author:
2004-02-11	2004-02-11	-	CMa

Q:

It sounds like the original problem (barcode scanner does not scan or searches too high for lower positions in towers) has re-surfaced. Perhaps there is an adjustment or some other option to repair this that Cosmas can suggest.

Is the 9 pin communication cable for the carousel a straight through cable or null-modem?

A:

It's a straight - other than the one for the handling. The scan height is in DM27. If you think that it does not scan at the proper height you may as well bend the BCR holder a little. Make sure that the scan line stays horizontal.