

Title:	Kiwis on same computer	ID:	
			131
17.07.13	Response:	Model:	Author:
		-	<i>CMa</i>

Kiwis on same computer

Test whether liconic STXGUI can be started twice with different databases to operate 2 Stt3k0 stores on the same workstation

on different windows xp computer (lab 3.2)

-> Firebird server installed

-> Firebird maestro installed, user ADMIN, pw test created

Liconic software copied to

D:\Data

a copy of the liconic software to

D:\Data2

D:\Data2\Liconic\Software\DriverConfig.ini file modified as follows:

[TCP]

port=**4333**

LogFile=D:\Data2\Liconic\STXTCPLLog.log

[WEB]

port=**82**

[XML]

port=**4334**

[paths]

StxMainFolder=D:\Data2\Liconic

CommunicationLogFolder=D:\Data2\Liconic\Log

LowLevelCommunicationLogFolder=D:\Data2\Liconic\LogLowLevel

XMLCommunicationLogFolder=D:\Data2\Liconic\XMLLog

PickJobReportFolder=D:\Data2\Liconic\PickJobReport

PickBackJobFolder=D:\Data2\Liconic\PickBackJob

[database]

server=localhost

path=D:\Data2\Liconic\DB\STXDRV.GDB

Embedded=0

[Options]

ShowErrorMessage=0

DevicesLogFile=D:\Data2\Liconic\DevError.log

KeepHistory=1

KeepEventsHistory=1

SimulationMode=0

PickingTubesVerify=1

STXGui.exe renamed to STXGui2.exe

→ Now StoreXGui could be started twice

