

Title:	How to Connect CO2	ID:	
		0130	
Date in:	Response:	Model:	Author:
2004-07-15	2004-07-15	STX200	CMa

Q:

A customer wants to run an incubator (Stx200, #3012) with Carbon Dioxide (it has not previously been run with CO2) by tomorrow afternoon. However- The incubator ICP has no tabs for running CO2. Could you let me know as soon as possible if this incubator model is fully enabled to run with CO2.

A:

This incubator was shipped with CO2.

Input pressure: ~2..4 bars.

Tube: 4x6mm

Read actual CO2 level: RD DM984 [1/100%]

Read CO2 level setting: RD DM894

Set CO2 level: WR DM894 xx (e.g. WR DM894 500 set CO2 to 5%)