

Title:	DC Unit Temperature too Low	ID:	
		0066	
Date in:	Response:	Model:	Author:
2006-04-27	2004-04-27	-	CMa

Q:

it looks as the fan never stops working on ser# 3318. Please prepare replacement parts for this unit, I will provide you with shipping address tomorrow


A:

On a DC the fan is connected to 513. It should either have the little add-on amplifier or - if there is an new Controller Card - go through the amplifier connected to the D-SUB 25Pin on the new Controller Card. There are two amplifiers on this card

Ampl1: Pin 1 (in) -> Pin 14 (out)

Ampl2: Pin 2 (in) -> Pin 15 (out)

The switched end of the fan is connected to Pin 14 of the 37Pin D-SUB of the Module (C1.3) (blue).



Add-On Ampl.

