

Title:	Castor STX SA Transfer Height	ID:	
		0028	
Date in:	Response:	Model:	Author:
2003-04-17	2003-04-17	StoreX SA Models	CMa

Q:

I notice from the photographs the STX40 is mounted on flat feet and the brochure indicates a pick off height of 1000mm.

Is there a revised height of 1030mm as per the STX1000 and 500 and is it also mounted on castors or are these possible options that must be specified at ordering

A:

The SA versions allow to lift the unit with a forklift. The lower subframe ~100mm is standard on SA's.

The transfer height is 1000mm but not very accurately. It can be adjusted with the feet.