

<b>Title:</b>	<b>Difference between STX BT and BTIT</b>			<b>ID:</b>
				0082
<b>Date in:</b>	<b>Response:</b>	<b>Model:</b>	<b>Author:</b>	
2006-03-17	2006-03-17	STX40	CMa	

**Q:**

We are working on the quotation for 4 incubators and we hesitate between STX 40 BT and STX 40 BT IT  
 On the list price, there is no description difference except the price  
 Could you explain what is the difference between the 2 models?  
 Thanks and best regards,

**A:**

The IT is mainly used for crystallization since the chiller is separated for the housing. The IT is also used when need reduced height. The standard unit is the STX40-ccBT.



