

Title:	Dimensions Incubation Chamber	ID:	
		<i>0084</i>	
Date in:	Response:	Model:	Author:
<i>2005-04-04</i>	<i>2005-04-04</i>	<i>STX500, 1000, STR602</i>	<i>CMa</i>

Q:

Can you provide the internal dimensions on the STX 500, 600 and 1000

A:

Dimensions WxDxH in mm

STX500: 615 x 600 x 1400

STX1000: 1120 x 640 x 1390

STX 620: 660 x 680 x 1620

Note that these dimensions can vary.

The width and height of the chamber is slightly larger (+~20mm) because of the design of the front of the chamber.