


Options

Easy Clean

Clean Incubation Chamber
Full Access for Cleaning
Simple Maintenance



Options

Easy Clean STX 40

With slide out stackers, the STX40 bottom surface is easily accessed for decontamination purpose. No other incubator has this option and it makes it impossible to clean their units under the robotic platform where the water is situated.

Order Information	
	<i>Order Nr.</i>
Easy Clean Option (STX40)	9118 12 21

Easy Clean STX220

In the event the incubator gets contaminated, the STX Series has the easy clean option that allows the carousel to be lifted into position to perform a full decontamination procedure.

No other incubator on the market has this option making it impossible to have the bottom of the unit cleaned.

Order Information	
	<i>Order Nr.</i>
Easy Clean Option (STX220)	9122 05 14